

Leading the way with new technologies and products

# MicroVention Product 2014 Catalog International



[microvention.com](http://microvention.com)

MicroVention, Inc.  
Worldwide Headquarters  
1311 Valencia Avenue  
Tustin, CA 92780 USA  
PH +1.714.247.8000  
PH 1.800.990.8368

MicroVention UK Limited  
Suite 3, The Barracks Building  
10 Cliffords Fort, North Shields  
Tyne and Wear, NE30 1JE  
United Kingdom  
PH +44 (0) 191 258 6777  
F +44 (0) 191 258 5999

MicroVention Europe, S.A.R.L.  
30 bis, rue du Vieil Abreuvoir  
78100 Saint-Germain-en-Laye  
France  
PH +33 (1) 39 21 77 46  
F +33 (1) 39 21 16 01

MicroVention Deutschland GmbH  
Hildebrandtstr. 4 F  
D-40215 Düsseldorf  
Germany  
PH +49 211 210 798-0  
F +49 211 210 798-29

# Table of Contents

## Products

### A. Coils

#### MicroPlex® Coil System

|                               |    |
|-------------------------------|----|
| - Cosmos® Complex Coils       | A1 |
| - Complex Platinum Coils      | A2 |
| - VFC® Complex Coils          | A3 |
| - HyperSoft® 3D Complex Coils | A4 |
| - Helical Platinum Coils      | A5 |
| - HyperSoft® Finishing Coils  | A6 |

#### Hydrogel Coated Coils

|  |    |
|--|----|
| - HydroFrame® HydroCoil Embolic System | A7 |
| - HydroFill® HydroCoil Embolic System  | A8 |
| - HydroSoft® HydroCoil Embolic System  | A9 |

### B. Delivery and Detachment Systems

|  |    |
|--|----|
| - V-Trak® Delivery Pusher and<br>V-Grip® Detachment Controller | B1 |
|--|----|

### C. Access Products

#### Guide Catheters

|                                     |    |
|-------------------------------------|----|
| - Chaperon® Guiding Catheter System | C1 |
| - Chaperon® Guiding Catheter Single | C2 |
| - SOFIA™ Distal Access Catheter     | C3 |

#### Guidewires

|  |    |
|--|----|
| - Traxcess® Guidewire and Docking Wire | C4 |
|--|----|

#### Microcatheters

|                              |    |
|------------------------------|----|
| - Headway® Microcatheter     | C5 |
| - Headway® Duo Microcatheter | C6 |

#### Balloon Catheters

|  |    |
|--|----|
| - Scepter C® Occlusion Balloon Catheter  | C7 |
| - Scepter XC® Occlusion Balloon Catheter | C8 |

### D. Intraluminal Devices

|   |    |
|---|----|
| - LVIS® Intraluminal Support Device     | D1 |
| - LVIS® Jr. Intraluminal Support Device | D2 |
| - FRED® System                          | D3 |

### E. Thrombectomy/Ischemic Stroke

|  |    |
|--|----|
| - ERIC® Retrieval Device and<br>ERIC® Retrieval Device Kit | E1 |
| - CASPER™ Carotid Artery Stent                             | E2 |

### F. Liquid Embolics

|                               |    |
|-------------------------------|----|
| - PHIL™ Liquid Embolic System | F1 |
|-------------------------------|----|

### G. Ordering Information

|               |    |
|---------------|----|
| International | G1 |
|---------------|----|

**Note:** this international catalog may contain products which are not available for sale or distribution in all countries, due to Regulatory requirements.



| <b>Cosmos® Complex Coils</b><br>For framing, filling and finishing                    |                            |                |
|---|----------------------------|----------------|
| 10 Coil Size  |                            |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen<br><b>(Stretch-Resistant)</b> |                            |                |
| V-Trak® Soft<br>Delivery System<br>Reference No.                                      | Large<br>Loop Dia.<br>(mm) | Length<br>(cm) |
| 100202CSSR-V*   | 2                          | 2              |
| 100254CSSR-V*   | 2.5                        | 4              |
| 100306CSSR-V*   | 3                          | 6              |
| 100408CSSR-V*   | 4                          | 8              |
| 100412CSSR-V*   | 4                          | 12             |
| 100515CSSR-V  | 5                          | 15             |
| 100522CSSR-V  | 5                          | 22             |
| 100618CSSR-V  | 6                          | 18             |
| 100626CSSR-V  | 6                          | 26             |
| 100722CSSR-V  | 7                          | 22             |
| 100731CSSR-V  | 7                          | 31             |
| 100825CSSR-V  | 8                          | 25             |
| 100837CSSR-V  | 8                          | 37             |
| 100933CSSR-V  | 9                          | 33             |
| 101036CSSR-V  | 10                         | 36             |

| <b>Cosmos® Complex Coils</b><br>For framing, filling and finishing |                            |                |
|--|----------------------------|----------------|
| 18 Coil Size   |                            |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen            |                            |                |
| V-Trak®<br>Delivery System<br>Reference No.                        | Large<br>Loop Dia.<br>(mm) | Length<br>(cm) |
| 180619CS-V   | 6                          | 19             |
| 180723CS-V   | 7                          | 23             |
| 180827CS-V   | 8                          | 27             |
| 180931CS-V   | 9                          | 31             |
| 181036CS-V   | 10                         | 36             |
| 181139CS-V   | 11                         | 39             |
| 181243CS-V   | 12                         | 43             |
| 181347CS-V   | 13                         | 47             |
| 181451CS-V   | 14                         | 51             |
| 181652CS-V   | 16                         | 52             |
| 181859CS-V   | 18                         | 59             |
| 182065CS-V   | 20                         | 65             |
| 182263CS-V   | 22                         | 63             |
| 182468CS-V   | 24                         | 68             |

\*2mm - 4mm features  Delivery System

## Features

Initial Progressively-Sized Loops

Combination of Unique  
Variable-Diameter Loops and  
Innovative Deployment

Stability

Choices of Length

## Benefits

Allows better controlled deployment into the aneurysm

Optimum framing provided by variable diameter loops directed towards periphery of aneurysm

As loops deploy, they brace the previously deployed loops

Loops find space along aneurysm periphery and provide excellent neck coverage

Stable basket upon subsequent coil deployment

Ability to concentrically fill

Long lengths provide greater volume fill and extensive coverage



| MicroPlex® Complex Coils                                |                      |                |
|---|----------------------|----------------|
| For framing and filling                                 |                      |                |
| 10 Coil Size  |                      |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen |                      |                |
| V-Trak®<br>Delivery System<br>Reference No.             | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
| 100204CC-V  | 2                    | 4              |
| 100307CC-V  | 3                    | 7              |
| 100410CC-V  | 4                    | 10             |
| 100512CC-V  | 5                    | 12             |
| 100615CC-V  | 6                    | 15             |
| 100718CC-V  | 7                    | 18             |
| 100820CC-V  | 8                    | 20             |
| 100924CC-V  | 9                    | 24             |
| 101026CC-V  | 10                   | 26             |

| MicroPlex® Complex Coils   |                      |                |
|--|----------------------|----------------|
| For framing and filling  |                      |                |
| 18 Coil Size   |                      |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen for coils 12mm or smaller<br>≥.018" Lumen for coils 13mm or larger |                      |                |
| 1 Diameter Loop  |                      |                |
| V-Trak®<br>Delivery System<br>Reference No.  | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
| 180410CC-V   | 4                    | 10             |
| 180512CC-V   | 5                    | 12             |
| 180615CC-V   | 6                    | 15             |
| 180717CC-V   | 7                    | 17             |
| 180820CC-V   | 8                    | 20             |
| 180923CC-V   | 9                    | 23             |
| 181026CC-V   | 10                   | 26             |
| 181128CC-V   | 11                   | 28             |
| 181231CC-V   | 12                   | 31             |
| 181332CC-V   | 13                   | 32             |
| 181434CC-V   | 14                   | 34             |
| 181639CC-V   | 16                   | 39             |
| 181844CC-V   | 18                   | 44             |
| 182050CC-V   | 20                   | 50             |

## Features

### Versatile Framework

## Benefits

Designed to reduce compartmentalization and provide greater neck coverage

Provides even distribution of coil loops at the periphery of the aneurysm

### Enhanced Stability

Each successive loop gently rotates course, changing direction as it is deployed

### Superior Conformability

Without assuming a pre-determined three-dimensional shape, the complex coil is able to frame a wide range of aneurysm morphologies

Facilitates placement of multiple consecutive complex coils



# Complex Coils



## VFC® Versatile Range Fill Coil

For filling

Recommended 2-Tip Marker Microcatheter  
≥ 0.0165" Lumen  
(Stretch-Resistant)

| V-Trak® Soft<br>Delivery System<br>Reference No. | Diameter<br>Range<br>(mm) | Length<br>(cm) |
|--|---------------------------|----------------|
| VFC010303-V*                                     | 1 - 3                     | 3              |
| VFC030606-V*                                     | 3 - 6                     | 6              |
| VFC030610-V*                                     | 3 - 6                     | 10             |
| VFC030615-V*                                     | 3 - 6                     | 15             |
| VFC061020-V                                      | 6 - 10                    | 20             |
| VFC061030-V                                      | 6 - 10                    | 30             |
| VFC101530-V                                      | 10 - 15                   | 30             |
| VFC101540-V                                      | 10 - 15                   | 40             |
| VFC152040-V                                      | 15 - 20                   | 40             |
| VFC152060-V                                      | 15 - 20                   | 60             |

\*1mm - 6mm features



## Features

Revolutionary Coil Design With  
Complex Loops And Waves

Single Coil Fills A Range  
Of Diameters

Long Length Options

Progressive Diameter Design

## Benefits

Allows for improved coil distribution, greater  
packing density and easier placement

Seeks out and fills open spaces

Efficient aneurysm coiling

Simplified coil selection

May significantly reduce coil inventory levels

Greater volume fill potential with a single coil

Less microcatheter movement

Softer coil during deployment

All coils are compatible with .0165" (0.42mm)  
ID microcatheters



### HyperSoft® 3D Complex Coil

For framing, filling and finishing

#### 10 Coil Size

Recommended 2-Tip Marker Microcatheter  
≥ 0.0165" Lumen

**Stretch-Resistant**

| V-Trak® Soft<br>Delivery System<br>Reference No. | Diameter<br>(mm) | Length<br>(cm) |
|--|------------------|----------------|
| 100102HS3D-V                                     | 1                | 2              |
| 100103HS3D-V                                     | 1                | 3              |
| 100152HS3D-V                                     | 1.5              | 2              |
| 100153HS3D-V                                     | 1.5              | 3              |
| 100154HS3D-V                                     | 1.5              | 4              |
| 100202HS3D-V                                     | 2                | 2              |
| 100203HS3D-V                                     | 2                | 3              |
| 100204HS3D-V                                     | 2                | 4              |
| 100254HS3D-V                                     | 2.5              | 4              |
| 100256HS3D-V                                     | 2.5              | 6              |
| 100304HS3D-V                                     | 3                | 4              |
| 100306HS3D-V                                     | 3                | 6              |
| 100355HS3D-V                                     | 3.5              | 5              |
| 100358HS3D-V                                     | 3.5              | 8              |
| 100406HS3D-V                                     | 4                | 6              |
| 100408HS3D-V                                     | 4                | 8              |
| 100510HS3D-V                                     | 5                | 10             |
| 100515HS3D-V                                     | 5                | 15             |

### Features

All-in-One Coil

Soft Complex Shape

Wider Range in Coil  
Sizes and Lengths

True to Size

### Benefits

Hybrid design provides ability to  
Frame, Fill, and Finish

Engineered to easily conform to  
various aneurysm morphologies

Provides optimal choices in complex  
aneurysm morphologies

Finish with Precision – Eliminates  
the need to under or oversize a coil

Featuring **V-Trak** Soft  
Delivery System



| MicroPlex® Helical Coils   |                      |                |
|--|----------------------|----------------|
| For filling  |                      |                |
| 10 Coil Size   |                      |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen<br>(Stretch-Resistant) |                      |                |
| V-Trak®<br>Delivery System<br>Reference No.                                    | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
| <b>Soft</b>  |                      |                |
| 100202HCSR-S-V   | 2                    | 2              |
| 100204HCSR-S-V   | 2                    | 4              |
| 100206HCSR-S-V   | 2                    | 6              |
| 100208HCSR-S-V   | 2                    | 8              |
| 100304HCSR-S-V   | 3                    | 4              |
| 100306HCSR-S-V   | 3                    | 6              |
| 100308HCSR-S-V   | 3                    | 8              |
| 100310HCSR-S-V   | 3                    | 10             |
| 100404HCSR-S-V   | 4                    | 4              |
| 100408HCSR-S-V   | 4                    | 8              |
| 100410HCSR-S-V   | 4                    | 10             |
| <b>Regular</b>   |                      |                |
| 100515HCSR-R-V   | 5                    | 15             |
| 100520HCSR-R-V   | 5                    | 20             |
| 100620HCSR-R-V   | 6                    | 20             |
| 100730HCSR-R-V   | 7                    | 30             |
| 100830HCSR-R-V   | 8                    | 30             |
| 100930HCSR-R-V   | 9                    | 30             |
| 101030HCSR-R-V   | 10                   | 30             |

| MicroPlex® Helical Coils                                |                      |                |
|---|----------------------|----------------|
| For filling   |                      |                |
| 18 Coil Size  |                      |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen |                      |                |
| V-Trak®<br>Delivery System<br>Reference No.             | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
| <b>Soft</b>   |                      |                |
| 180204HC-S-V  | 2                    | 4              |
| 180208HC-S-V  | 2                    | 8              |
| 180304HC-S-V  | 3                    | 4              |
| 180308HC-S-V  | 3                    | 8              |
| 180406HC-S-V  | 4                    | 6              |
| 180410HC-S-V  | 4                    | 10             |
| 180512HC-S-V  | 5                    | 12             |
| 180615HC-S-V  | 6                    | 15             |
| <b>Regular</b>  |                      |                |
| 180520HC-R-V  | 5                    | 20             |
| 180620HC-R-V  | 6                    | 20             |
| 180730HC-R-V  | 7                    | 30             |
| 180830HC-R-V  | 8                    | 30             |
| 180930HC-R-V  | 9                    | 30             |
| 181030HC-R-V  | 10                   | 30             |
| 181230HC-R-V  | 12                   | 30             |
| 181430HC-R-V  | 14                   | 30             |
| 181630HC-R-V  | 16                   | 30             |
| 181830HC-R-V  | 18                   | 30             |
| 182030HC-R-V  | 20                   | 30             |

## Features

Versatile Framework

Available in Wide Range of Sizes

## Benefits

Exceptionally soft for easy deployment and repositioning

10 and 18 systems with 10 system stretch-resistant available

Provides efficient packing



## MicroPlex® HyperSoft® Coils

For finishing

### 10 Coil Size

Recommended 2-Tip Marker Microcatheter  
≥.0165" Lumen  
(Stretch-Resistant)

| V-Trak® Soft<br>Delivery System<br>Reference No. | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
|--|----------------------|----------------|
| 100151HS-V                                       | 1.5                  | 1              |
| 100152HS-V                                       | 1.5                  | 2              |
| 100153HS-V                                       | 1.5                  | 3              |
| 100154HS-V                                       | 1.5                  | 4              |
| 100201HS-V                                       | 2                    | 1              |
| 100202HS-V                                       | 2                    | 2              |
| 100203HS-V                                       | 2                    | 3              |
| 100204HS-V                                       | 2                    | 4              |
| 100206HS-V                                       | 2                    | 6              |
| 100254HS-V                                       | 2.5                  | 4              |
| 100256HS-V                                       | 2.5                  | 6              |
| 100304HS-V                                       | 3                    | 4              |
| 100306HS-V                                       | 3                    | 6              |
| 100308HS-V                                       | 3                    | 8              |
| 100404HS-V                                       | 4                    | 4              |
| 100406HS-V                                       | 4                    | 6              |
| 100408HS-V                                       | 4                    | 8              |
| 100506HS-V                                       | 5                    | 6              |
| 100508HS-V                                       | 5                    | 8              |
| 100606HS-V                                       | 6                    | 6              |
| 100608HS-V                                       | 6                    | 8              |

## Features

Soft and Conformable Coil

Softer Pusher

New 1.5mm and 2.5mm Sizes

## Benefits

Greater packing density and easier placement

Ability to place coil in confined spaces

Less microcatheter movement or deflection for better control placing the coil into the aneurysm

Designed to reduce microcatheter kickback

1.5mm with 4 lengths and 2.5 with 2 lengths offer greater flexibility of choice

Featuring **V-Trak** Delivery System  
*Soft*





| <b>HydroFrame® Complex Coils</b><br>For framing and filling                    |                      |                |
|--|----------------------|----------------|
| <b>10 Coil Size</b>  |                      |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen<br>(Stretch-Resistant) |                      |                |
| V-Trak® Soft<br>Delivery System<br>Reference No.                               | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
| 100202HFRM-V*  | 2                    | 2              |
| 100254HFRM-V*  | 2.5                  | 4              |
| 100306HFRM-V*  | 3                    | 6              |
| 100405HFRM-V*  | 4                    | 5              |
| 100408HFRM-V*  | 4                    | 8              |
| 100510HFRM-V   | 5                    | 10             |
| 100515HFRM-V   | 5                    | 15             |
| 100612HFRM-V   | 6                    | 12             |
| 100619HFRM-V   | 6                    | 19             |
| 100715HFRM-V   | 7                    | 15             |
| 100728HFRM-V   | 7                    | 28             |
| 100817HFRM-V   | 8                    | 17             |
| 100833HFRM-V   | 8                    | 33             |
| 100931HFRM-V   | 9                    | 31             |
| 101036HFRM-V   | 10                   | 36             |

| <b>HydroFrame® Complex Coils</b><br>For framing and filling |                   |                |
|---|-------------------|----------------|
| <b>18 Coil Size</b>   |                   |                |
| Recommended 2-Tip Marker Microcatheter<br>≥.0165" Lumen     |                   |                |
| V-Trak®<br>Delivery System<br>Reference No.                 | Loop Dia.<br>(mm) | Length<br>(cm) |
| 180619HFRM-V  | 6                 | 19             |
| 180723HFRM-V  | 7                 | 23             |
| 180827HFRM-V  | 8                 | 27             |
| 180931HFRM-V  | 9                 | 31             |
| 181036HFRM-V  | 10                | 36             |
| 181139HFRM-V  | 11                | 39             |
| 181243HFRM-V  | 12                | 43             |
| 181347HFRM-V  | 13                | 47             |
| 181445HFRM-V  | 14                | 45             |
| 181644HFRM-V  | 16                | 44             |
| 181850HFRM-V  | 18                | 50             |
| 182048HFRM-V  | 20                | 48             |

\*2mm - 4mm features **V-Trak® Soft** Delivery System

## Features

First Two Loops Progressively-Sized

Variable-Diameter Loops and  
Innovative Deployment

Stability

Ease of Use

Greater Volume Fill

Ability to Deliver Hydrogel  
at the Neck

## Benefits

Better initial deployment

Excellent framing along aneurysm periphery and  
neck

Stable basket upon subsequent coil deployment

Ability to concentrically fill

No prepping or steaming

HydroFrame provides greater volume fill than other  
10 system framing coils

Hydrogel at the neck has been shown to reduce  
recurrence rates<sup>1</sup>

1. H.J. Cloft, for the HEAL Investigators. HydroCoil for Endovascular Aneurysm Occlusion (HEAL) Study: 3-6 Month Angiographic Follow-Up Results. AJNR AM J Neuroradiology; 28:152-154, Jan 2007 (197 aneurysms; 135 follow-up cases)



## HydroFill® Helical Coils

For filling

Recommended 2-Tip Marker Microcatheter  
≥.0165" Lumen  
(Stretch-Resistant)

| V-Trak® Soft<br>Delivery System<br>Reference No. | Loop<br>Dia.<br>(mm) | Length<br>(cm) |
|--|----------------------|----------------|
| 100202HFIL-V*                                    | 2                    | 2              |
| 100204HFIL-V*                                    | 2                    | 4              |
| 100206HFIL-V*                                    | 2                    | 6              |
| 100304HFIL-V*                                    | 3                    | 4              |
| 100306HFIL-V*                                    | 3                    | 6              |
| 100408HFIL-V*                                    | 4                    | 8              |
| 100410HFIL-V*                                    | 4                    | 10             |
| 100510HFIL-V                                     | 5                    | 10             |
| 100515HFIL-V                                     | 5                    | 15             |
| 100610HFIL-V                                     | 6                    | 10             |
| 100615HFIL-V                                     | 6                    | 15             |
| 100620HFIL-V                                     | 6                    | 20             |
| 100715HFIL-V                                     | 7                    | 15             |
| 100720HFIL-V                                     | 7                    | 20             |
| 100820HFIL-V                                     | 8                    | 20             |
| 100830HFIL-V                                     | 8                    | 30             |
| 100920HFIL-V                                     | 9                    | 20             |
| 100930HFIL-V                                     | 9                    | 30             |
| 101020HFIL-V                                     | 10                   | 20             |
| 101030HFIL-V                                     | 10                   | 30             |
| 101230HFIL-V                                     | 12                   | 30             |
| 101440HFIL-V                                     | 14                   | 40             |
| 101640HFIL-V                                     | 16                   | 40             |
| 101850HFIL-V                                     | 18                   | 50             |
| 102050HFIL-V                                     | 20                   | 50             |
| 102250HFIL-V                                     | 22                   | 50             |
| 102450HFIL-V                                     | 24                   | 50             |

## Features

No Steaming or Preparation  
of Coil Required

Longer Working Time

Volume Fill

Longer Lengths

Hydrogel-Based Coil System

Ability to Perform All-Hydrogel  
Cases

V-Trak® Soft Delivery System,  
Sizes 2mm - 4mm

## Benefits

Saves time and work of preparing coil

Ability to utilize more time if needed  
to deploy coil:

**Microcatheters ≥ .0165"**

- 30 minutes for 2mm - 4mm
- 10 minutes for 5mm - 24mm

**Microcatheters ≥ .021"**

- 30 minutes for all sizes

2mm - 4mm .016" expanded diameter  
5mm - 24mm .018" expanded diameter

Up to 50cm on larger sizes

Clinical benefits of hydrogel, proven in  
numerous clinical studies, including HELPS

Ability to utilize easy-to use hydrogel  
coils for framing, filling and finishing:  
HydroFrame®, HydroFill®, and  
HydroSoft® Embolic Coil Systems

Less microcatheter movement or deflection  
for better control placing the coil into the  
aneurysm

Designed to reduce microcatheter kickback

\*2mm - 4mm features





## HydroSoft® Helical Coils

For filling and finishing

### 10 Coil Size

Recommended 2-Tip Marker Microcatheter  
≥.0165" Lumen  
(Stretch-Resistant)

| V-Trak® Soft<br>Delivery System<br>Reference No. | Loop Dia.<br>(mm) | Length<br>(cm) |
|--|-------------------|----------------|
| 100152H2HS-V*                                    | 1.5               | 2              |
| 100154H2HS-V*                                    | 1.5               | 4              |
| 100201H2HS-V*                                    | 2                 | 1              |
| 100202H2HS-V*                                    | 2                 | 2              |
| 100203H2HS-V*                                    | 2                 | 3              |
| 100204H2HS-V*                                    | 2                 | 4              |
| 100206H2HS-V*                                    | 2                 | 6              |
| 100208H2HS-V*                                    | 2                 | 8              |
| 100254H2HS-V*                                    | 2.5               | 4              |
| 100256H2HS-V*                                    | 2.5               | 6              |
| 100304H2HS-V*                                    | 3                 | 4              |
| 100306H2HS-V*                                    | 3                 | 6              |
| 100308H2HS-V*                                    | 3                 | 8              |
| 100310H2HS-V*                                    | 3                 | 10             |
| 100404H2HS-V*                                    | 4                 | 4              |
| 100406H2HS-V*                                    | 4                 | 6              |
| 100408H2HS-V*                                    | 4                 | 8              |
| 100410H2HS-V*                                    | 4                 | 10             |
| 100506H2HS-V                                     | 5                 | 6              |
| 100508H2HS-V                                     | 5                 | 8              |
| 100510H2HS-V                                     | 5                 | 10             |
| 100515H2HS-V                                     | 5                 | 15             |
| 100520H2HS-V                                     | 5                 | 20             |
| 100606H2HS-V                                     | 6                 | 6              |
| 100608H2HS-V                                     | 6                 | 8              |
| 100610H2HS-V                                     | 6                 | 10             |
| 100615H2HS-V                                     | 6                 | 15             |

## Features

Gel on the Inside

Microcatheter Compatibility

Stretch-Resistant

## Benefits

No prep, no steaming

Provides the softness of a platinum finishing coil

30 minute working time and compatible with .0165" - .021" microcatheters

Facilitates repositioning

\*1.5mm - 4mm features



Delivery System

### V-Trak Delivery Pusher

#### Features

No Preparation

More Flexible Delivery Pusher

Softer, Shorter  
Detachment Zone

Pusher Stability

Coil Exit Markers

#### Benefits

Simplifies operation; reduces procedure time

Minimizes microcatheter movement;  
enhances coil deliverability

Superior coil placement ability;  
reduces microcatheter tip deflection

Solid mandrel pusher offers more  
kink-resistance and improved coil  
advancement

Indicates coil is approaching  
microcatheter tip, reducing  
fluoroscopy time



### V-Grip

#### Features

Self-contained Integrated  
Power Supply

Rapid Coil Detachment

Sterile, Disposable Design

Ergonomic Design

Visual and Audible Indicators

#### Benefits

Ready to use; no accessories, no patient  
grounding, no connecting cables or  
additional batteries

Single person operation; simple one-button  
coil detachment in 0.75 seconds for maximum  
efficiency and speed

Reduces set-up and procedure time

Lightweight, no device preparation, and  
easily identifiable

Provides instant notification that  
device is ready for use

#### V-Grip<sup>®</sup> Detachment Controller

For use with the V-Trak<sup>®</sup> Delivery System

**5 Detachment Controllers Per Box**

**Reference No.**

VG501



### Chaperon® Guiding Catheter System

Guiding Catheter System

| 5F Neuro Guiding Catheter / System |        |                         |         |             | 6F Neuro Guiding Catheter / System |        |                         |         |             |
|------------------------------------|--------|-------------------------|---------|-------------|------------------------------------|--------|-------------------------|---------|-------------|
| Product Code                       | Type   | OD Inner/Outer (French) | Shape   | ID (inches) | Product Code                       | Type   | OD Inner/Outer (French) | Shape   | ID (inches) |
| GC595STJB                          | System | 4/5                     | STR/JB2 | 0.059       | GC695STJB                          | System | 5/6                     | STR/JB2 | 0.071       |
| GC595STVT                          | System | 4/5                     | STR/VTR | 0.059       | GC695STVT                          | System | 5/6                     | STR/VTR | 0.071       |
| GC595STSI                          | System | 4/5                     | STR/SIM | 0.059       | GC695STSI                          | System | 5/6                     | STR/SIM | 0.071       |
| GC595M2JB                          | System | 4/5                     | MP2/JB2 | 0.059       | GC695M2JB                          | System | 5/6                     | MP2/JB2 | 0.071       |
| GC595M2SI                          | System | 4/5                     | MP2/SIM | 0.059       | GC695M2SI                          | System | 5/6                     | MP2/SIM | 0.071       |
| GC595M2VT                          | System | 4/5                     | MP2/VTR | 0.059       | GC695M2VT                          | System | 5/6                     | MP2/VTR | 0.071       |
|                                    |        |                         |         |             | GC695BUJB                          | System | 5/6                     | BUR/JB2 | 0.071       |
|                                    |        |                         |         |             | GC695BUVT                          | System | 5/6                     | BUR/VTR | 0.071       |

#### Guiding Catheter Shapes



STR



MP2



BUR

#### Inner Catheter Shapes



VTR



JB2



SIM2

## Features

Inner Catheter Combined System

Largest ID

Flat Braid Reinforced Shaft


## Benefits

Atraumatic access and deeper positioning for more support with high maximum injection pressure (**4F**-750 PSI **5F**-1000 PSI)

Enhanced work space for multiple catheter usage (**5F**-0.059" **6F**-0.071")

Increased kink resistance and greater torque control

| <b>Chaperon® Guiding Catheter Single</b><br>Guiding Catheter / Single |        |                   |       |             |   |        |                   |       |             |
|---|--------|-------------------|-------|-------------|---|--------|-------------------|-------|-------------|
| <b>5F</b> Neuro Guiding Catheter / Single                             |        |                   |       |             | <b>6F</b> Neuro Guiding Catheter / Single |        |                   |       |             |
| Product Code  | Type   | OD Outer (French) | Shape | ID (inches) | Product Code                              | Type   | OD Outer (French) | Shape | ID (inches) |
| <b>GC595ST</b>  | Single | 5                 | STR   | 0.059       | <b>GC695ST</b>                            | Single | 6                 | STR   | 0.071       |
| <b>GC595M2</b>  | Single | 5                 | MP2   | 0.059       | <b>GC695M2</b>                            | Single | 6                 | MP2   | 0.071       |



## Features

Largest ID

Flat Braid Reinforced Shaft

Proprietary Polyester for Outer Material

## Benefits

Enhanced work space for multiple catheter usage  
(**5F**-0.059" **6F**-0.071")

Increased kink resistance and greater torque control

Provides better support

| <b>SOFIA™ Distal Access Catheter</b> Soft Torqueable Catheter Optimized For Intracranial Access |                        |                       |                     |           |                     |                    |                  |
|---|------------------------|-----------------------|---------------------|-----------|---------------------|--------------------|------------------|
| 1 per box and includes shaping mandrel and introducer sheath                                    |                        |                       |                     |           |                     |                    |                  |
| Product Code  | Catheter Size (French) | Proximal OD (inch/mm) | Distal OD (inch/mm) | ID (inch) | Working Length (cm) | Distal Length (cm) | Distal Tip Shape |
| DA5125ST  | 5F                     | 0.068 / 1.7           | 0.068 / 1.7         | 0.055     | 125                 | 17                 | STRAIGHT         |

## Features

Exceptionally Soft Distal Shaft

Low Profile

Steam Shapeable Tip and Torqueable Shaft

Enhanced Kink Resistance

## Benefits

Easier navigation in tortuous vessels

Compatible with 6F .070 inch ID guide catheter systems

Steerable around bifurcations

Easy to handle with 1:1 push/pull control



# Guidewire with Hydrophilic Coating

| Traxcess® Guidewire and Docking Wire   |              |             |           |                          |                             |                           |
|--|--------------|-------------|-----------|--------------------------|-----------------------------|---------------------------|
| Guidewire (0.014")   |              |             |           |                          |                             |                           |
| One unit per box / Guidewires include shaping mandrel, insertion tool, and torque device |              |             |           |                          |                             |                           |
| Product Name   | Product Code | Description | Tip Shape | Total/Distal Length (cm) | Proximal/Distal OD (inches) | Shapeable Tip Length (cm) |
| Traxcess 14  | GW1420040    | Soft        | STR       | 200 / 40                 | 0.014 / 0.012               | 1.4                       |
| Traxcess 14EX  | GW1420040X   | Support     | STR       | 200 / 40                 | 0.014 / 0.012               | 1.4                       |
| Docking Wire (0.014")  |              |             |           |                          |                             |                           |
| Traxcess Docking Wire  | GW14100EX    | Docking     | —         | 115                      | 0.014                       | —                         |



## Features

- Nitinol Property Control Technology
- Hybrid Guidewire Core Construction
- Terumo Glide Technology™

## Benefits

- Excellent tip shapability and shape retention
- 1:1 torque transmission and direct wire control
- Enhanced durability and lubricity in hydrophilic coating



## Features

- Compatible with Traxcess Guidewires
- Simple Docking Wire System

## Benefits

- Total length is 313cm when combined with the Traxcess Guidewire (2cm overlap)
- Easily convert the Traxcess Guidewire into exchange length



## Headway® Microcatheter

Microcatheter with Hydrophilic Coating

One unit per box and includes shaping mandrel and introducer sheath

| Product Name              | Description       | Product Code | ID (inches) | Tip Shape | Usable Length (cm) | OD Prox./Distal (French) | Tip Markers |
|---------------------------|-------------------|--------------|-------------|-----------|--------------------|--------------------------|-------------|
| Headway17 Advanced 'Soft' | Straight          | MC172150S    | 0.017       | STR       | 150                | 2.4 / 1.7                | 2           |
| Headway17 Advanced        | Straight          | MC172150SX   | 0.017       | STR       | 150                | 2.4 / 1.7                | 2           |
| Headway17 Advanced        | Pre-shaped tip 45 | MC17215045X  | 0.017       | 45        | 150                | 2.4 / 1.7                | 2           |
| Headway17 Advanced        | Pre-shaped tip 90 | MC17215090X  | 0.017       | 90        | 150                | 2.4 / 1.7                | 2           |
| Headway17 Advanced        | Pre-shaped tip J  | MC172150JX   | 0.017       | J         | 150                | 2.4 / 1.7                | 2           |
| Headway21                 | Straight          | MC212150S    | 0.021       | STR       | 150                | 2.5 / 2.0                | 2           |
| Headway27                 | Straight          | MC272150S    | 0.027       | STR       | 150                | 3.0 / 2.6                | 2           |

## Product

### Headway17 Advanced

0.017 inch ID microcatheter with atraumatic 1.7F distal tip with enhanced catheter control, superior shapability and shape retention. An enhanced distal shaft provides excellent support for coiling while the segmented design ensures smooth trackability

### Headway17 Advanced 'Soft'

0.017 inch ID microcatheter with a softer distal shaft to provide more sensitivity to coil deployment force

### Headway21

0.021 inch ID microcatheter with a 2.0F rounded distal tip, optimized for stent delivery and the ability to track through tortuous anatomy

### Headway27

0.027 inch ID microcatheter with 2.6F distal tip diameter, designed specifically for flow diverting stent delivery

## Benefits



# Innovative Catheter Technology

## Headway® Duo Microcatheter

Microcatheter with Hydrophilic Coating

1 per box; includes shaping mandrel and introducer sheath

| Product Name       | Product Code | ID Body/Tip (inches) | Tip Shape | Usable Length (cm) | Distal Length (cm) | OD Prox./Dist. (French) | Tip Markers |
|--------------------|--------------|----------------------|-----------|--------------------|--------------------|-------------------------|-------------|
| Headway Duo, 156cm | MC162156S    | 0.0165 / 0.0165      | STR       | 156                | 9                  | 2.1 / 1.6               | 2           |
| Headway Duo, 167cm | MC162167S    | 0.0165 / 0.013       | STR       | 167                | 9                  | 2.1 / 1.3               | 1           |

## Features

### Hybrid Braid/Coil Design

Tight pitch coil provides lumen integrity and bending flexibility  
Proximal braid provides firmness and a torqueable catheter body  
Will not ovalize or kink

### Small OD, Large ID

0.0165" body ID microcatheter with low profile distal tip  
156cm Headway® Duo microcatheter compatible with coils and  $\leq 0.014$ " guidewires  
167cm Headway® Duo microcatheter compatible with  $\leq 0.014$ " /  $0.012$ " (proximal/distal) guidewires

### Superior Trackability

True 1:1 push/pull while tracking through tortuous anatomy

### Liquid Embolic Compatibility

Compatible with Onyx™ Liquid Embolic and nBCA  
High burst strength  
Long length allows distal lesion access

### Soft Distal Tip

Smooth, atraumatic tracking  
High sensitivity to coil deployment force  
Excellent tip "painting" ability

## Scepter C® Occlusion Balloon Catheter

Compliant Balloon Catheter



| Product Code               | Description                       | Balloon Diameter (mm) | Balloon Length (mm) | Distal Tip Length (mm) |
|----------------------------|-----------------------------------|-----------------------|---------------------|------------------------|
| 1 balloon catheter per box |                                   |                       |                     |                        |
| BC0410C                    | Compliant Occlusion Balloon 4.0mm | 4.0                   | 10                  | 5                      |
| BC0415C                    | Compliant Occlusion Balloon 4.0mm | 4.0                   | 15                  | 5                      |
| BC0420C                    | Compliant Occlusion Balloon 4.0mm | 4.0                   | 20                  | 5                      |

### Features

#### 14 Wire Compatible

Provides more choices of guidewire (Traxcess® docking wire compatible)  
Improves trackability and stability with better performance of the wire

#### Hydrophilic Coating on the Balloon

Significantly improves trackability  
Enables use of longer balloons without compromising trackability

#### Soft & Long Tip

Improves balloon stability and trackability

#### Low Profile

Ensures compatibility with double catheter technique in 6F guide  
Improves trackability

#### Distal Tip Marker

Positive identification of distal tip location

## Scepter<sup>XC</sup>® Occlusion Balloon Catheter

X-tra Compliant Balloon Catheter



| Product Code               | Description                                    | Balloon Diameter (mm) | Balloon Length (mm) | Distal Tip Length (mm) |
|----------------------------|--|-----------------------|---------------------|------------------------|
| 1 balloon catheter per box |  |                       |                     |                        |
| BC0411XC                   | X-tra Compliant Occlusion Balloon 4.0mm x 11mm | 4.0                   | 11                  | 5                      |

### Features

#### 14 Wire Compatible

Provides more choices of guidewire (Traxcess® docking wire compatible)  
Improves trackability and stability with better performance of the wire

#### Hydrophilic Coating on the Balloon

Significantly improves trackability  
Enables use of longer balloons without compromising trackability

#### Soft & Long Tip

Improves balloon stability and trackability

#### Low Profile

Ensures compatibility with double catheter technique in 6F guide  
Improves trackability

#### Distal Tip Marker

Positive identification of distal tip location

## LVIS® Device

Intraluminal Support Device

Low-profile Visualized Intraluminal Support

One unit per box

| Product Code | Unconstrained Diameter x Total Length (mm) | Total Length*/Working Length in Different Vessel Diameters (mm) |         |         |         |         |         |         |
|--------------|--|---|---------|---------|---------|---------|---------|---------|
|              |  | 2.5mm   | 3.0mm   | 3.5mm   | 4.0mm   | 4.5mm   | 5.0mm   | 5.5mm   |
| 212517-CAS   | 3.5 x 19                                   | 22 / 18   | 20 / 16 | 19 / 15 |         |         |         |         |
| 212525-CAS   | 3.5 x 24                                   | 30 / 26   | 27 / 23 | 24 / 20 |         |         |         |         |
| 213015-CAS   | 4.5 x 16                                   |   | 20 / 16 | 19 / 15 | 18 / 14 | 16 / 12 |         |         |
| 213025-CAS   | 4.5 x 23                                   |   | 30 / 26 | 28 / 24 | 26 / 22 | 23 / 19 |         |         |
| 213041-CAS   | 4.5 x 34                                   |   | 46 / 42 | 43 / 39 | 39 / 35 | 34 / 30 |         |         |
| 214035-CAS   | 5.5 x 27                                   |   |         |         | 40 / 36 | 37 / 33 | 32 / 28 | 27 / 23 |
| 214049-CAS   | 5.5 x 34                                   |   |         |         | 54 / 50 | 49 / 45 | 43 / 39 | 34 / 30 |

Notes: Compatible Microcatheter: Headway®21

\*Total Length (includes flared ends) = Working Length + 4mm

## Features

Woven braid device made with 16 passes of a single nitinol wire

1mm compliant cell design

Radiopaque markers at both ends and two radiopaque helical strands along the length of the stent

Up to 75% retractable

## Benefits

Up to 22% metal coverage provides effective scaffold across aneurysm neck and optimal vessel wall apposition

Allows for easier passage of microcatheters

Enhanced radiopacity helps ensure optimal deployment and visibility

Improves accuracy in targeting the landing zone during deployment



# Intraluminal Support Device

## LVIS® Jr. Device

Intraluminal Support Device

Low-profile Visualized Intraluminal Support

One unit per box

| Product Code | Unconstrained Diameter<br>x<br>Total Length<br>(mm) | Total Length* / Working Length<br>in Different Vessel Diameters<br>(mm) |         |         |         |
|--------------|---|---|---------|---------|---------|
|              |   | 2.0mm   | 2.5mm   | 3.0mm   | 3.5mm   |
| 172010-CASJ  | 2.5 x 13  | 14 / 10   | 13 / 9  |         |         |
| 172014-CASJ  | 2.5 x 17  | 18 / 14   | 17 / 13 |         |         |
| 172020-CASJ  | 2.5 x 23  | 24 / 20   | 23 / 19 |         |         |
| 172032-CASJ  | 2.5 x 34  | 36 / 32   | 34 / 30 |         |         |
| 172516-CASJ  | 3.5 x 18  |   | 20 / 16 | 19 / 15 | 18 / 14 |
| 172524-CASJ  | 3.5 x 23  |   | 27 / 23 | 25 / 21 | 23 / 19 |
| 172530-CASJ  | 3.5 x 28  |   | 34 / 30 | 32 / 28 | 28 / 24 |
| 172537-CASJ  | 3.5 x 33  |   | 40 / 36 | 37 / 33 | 33 / 29 |

Notes: Compatible Microcatheter: Headway®17 or Scepter Occlusion Balloon

\*Total Length (includes flared ends) = Working Length + 4mm

## Features

Compatible with Headway®17 Advanced microcatheter

Deliverable through Scepter Occlusion Balloon Catheter

1.5mm compliant cell design

Up to 75% retractable

## Benefits

Low profile 1.7F delivery system improves navigability

Post-balloon stenting provides new method of stent utilization

Allows for easier passage of microcatheters

Improves accuracy in targeting the landing zone during deployment

| <b>FRED® System</b><br>Flow Re-Direction Endoluminal Device |                  |                                   |                                      |   |
|---|------------------|-----------------------------------|--------------------------------------|---|
| One unit per box  |                  |                                   |                                      |   |
| Cross-Reference Product Code                                | New Product Code | Fully Open Implant Diameter* (mm) | Total Length / Working Length** (mm) | Fully Open Total Length / Working Length** (mm) |
|   |                  |                                   | 3.0mm OD                             | 3.5mm OD  |
| FRED3505  | FRED3507         | 3.5                               | 16 / 10                              | 13 / 7  |
| N/A   | FRED3511         | 3.5                               | 22 / 17                              | 17 / 11   |
| FRED3525  | FRED3516         | 3.5                               | 30 / 24                              | 22 / 16   |
| N/A   | FRED3524         | 3.5                               | 40 / 34                              | 31 / 24   |
| N/A   | FRED3536         | 3.5                               | 54 / 49                              | 40 / 36   |
|   |                  |                                   | 3.5mm OD                             | 4.0mm OD  |
| N/A   | FRED4007         | 4.0                               | 16 / 10                              | 13 / 7  |
| FRED4015  | FRED4012         | 4.0                               | 24 / 18                              | 18 / 12   |
| FRED4025  | FRED4017         | 4.0                               | 31 / 26                              | 23 / 17   |
| N/A   | FRED4026         | 4.0                               | 43 / 37                              | 32 / 26   |
| N/A   | FRED4038         | 4.0                               | 58 / 53                              | 44 / 38   |
|   |                  |                                   | 4.0mm OD                             | 4.5mm OD  |
| N/A   | FRED4508         | 4.5                               | 18 / 11                              | 15 / 8  |
| FRED4515  | FRED4513         | 4.5                               | 26 / 19                              | 20 / 13   |
| FRED4525  | FRED4518         | 4.5                               | 34 / 27                              | 25 / 18   |
| N/A   | FRED4528         | 4.5                               | 45 / 39                              | 34 / 28   |
| N/A   | FRED4539         | 4.5                               | 62 / 56                              | 45 / 39   |
|   |                  |                                   | 4.5mm OD                             | 5.0mm OD  |
| N/A   | FRED5009         | 5.0                               | 18 / 12                              | 15 / 9  |
| FRED5015  | FRED5014         | 5.0                               | 26 / 20                              | 21 / 14   |
| FRED5025  | FRED5019         | 5.0                               | 35 / 28                              | 26 / 19   |
| N/A   | FRED5029         | 5.0                               | 48 / 41                              | 36 / 29   |
|   |                  |                                   | 5.0mm OD                             | 5.5mm OD  |
| FRED5515  | FRED5514         | 5.5                               | 28 / 22                              | 22 / 14   |
| FRED5535  | FRED5526         | 5.5                               | 46 / 38                              | 32 / 26   |

## Features

### Integrated Dual-Layer Coverage

## Benefits

16-wire outer stent with 48-wire inner stent for enhanced radial force  
Inner stent provides high percentage metal surface area for flow-redirection

### Braided Self-Expanding Nickel Titanium, Paired-Stent Design

Enables simultaneous, dynamic push/pull dual-stent deployment with single operator ease of use  
Easy stent opening and superb vessel apposition, with partial resheathability

### Radiopaque Distal Wire Tip, Proximal and Distal Markers, and Interweaved Working Length Helical Strands


Full-stent length fluoroscopic visualization

Notes:

Compatible with Headway® 27 Microcatheter

\* Max OD constrained in vessel

\*\* Outer device total length (including flared ends) and inner device working length (marked by helical strands), constrained in vessel





# Retrieval Device

## ERIC® Retrieval Device

Embolus Retriever with Interlinked Cages

| Product Code                             | Product Name | # of Spheres | Diameter (mm) | Working Length (mm) |
|--|--------------|--------------|---------------|---------------------|
| 1 per box and includes introducer sheath |              |              |               |                     |
| ER173015                                 | ERIC 3       | 3            | 3             | 15                  |
| ER173020                                 | ERIC 3       | 4            | 3             | 20                  |
| ER174024                                 | ERIC 4       | 4            | 4             | 24                  |
| ER174030                                 | ERIC 4       | 5            | 4             | 30                  |
| ER176044                                 | ERIC 6       | 5            | 6             | 44                  |

## Features

### Interlinked Cage Design

## Benefits

- Allows immediate clot-device integration at deployment
- Independent cages that conform to tight curves to retain clot
- Minimum contact to the vessel wall for atraumatic retrieval

### Headway17 Advanced Compatible

- Improves access in the tortuous anatomy
- Leaves more space for more suction with a small outer profile
- Working length can be adjusted by re-sheathing

## ERIC® Retrieval Device Kit

Embolus Retriever with Interlinked Cages

Comes with  
**HeadwayAdvanced**  
Microcatheter

| Product Code  | Product Name | Description                        |
|---|--------------|------------------------------------|
| Includes 1 ERIC Retrieval Device and 1 Headway17 Advanced Microcatheter |              |                                    |
| ER315KTA  | ERIC 3 Kit   | ERIC 3x15 with Headway 17 Advanced |
| ER320KTA  | ERIC 3 Kit   | ERIC 3x20 with Headway 17 Advanced |
| ER424KTA  | ERIC 4 Kit   | ERIC 4x24 with Headway 17 Advanced |
| ER430KTA  | ERIC 4 Kit   | ERIC 4x30 with Headway 17 Advanced |
| ER644KTA  | ERIC 6 Kit   | ERIC 6x44 with Headway 17 Advanced |



## CASPER<sup>TM</sup> RX Carotid Artery Stent

1 per box

| Product Code   | Unconstrained Dimensions |                        |                     | Implanted Dimensions         |                     |                              |                     |
|----------------|--------------------------|------------------------|---------------------|------------------------------|---------------------|------------------------------|---------------------|
|                | Diameter (mm)            | Dual Layer Length (mm) | Overall Length (mm) | 1mm < Unconstrained Diameter |                     | 2mm < Unconstrained Diameter |                     |
|                |                          |                        |                     | Dual Layer Length (mm)       | Overall Length (mm) | Dual Layer Length (mm)       | Overall Length (mm) |
| CPR-0520-143RX | 5                        | 20                     | 25                  | 20                           | 33                  | 22                           | 35                  |
| CPR-0530-143RX | 5                        | 30                     | 37                  | 35                           | 47                  | 38                           | 52                  |
| CPR-0540-143RX | 5                        | 40                     | 47                  | 45                           | 59                  | 52                           | 64                  |
| CPR-0616-143RX | 6                        | 16                     | 22                  | 20                           | 32                  | 23                           | 35                  |
| CPR-0625-143RX | 6                        | 25                     | 33                  | 30                           | 44                  | 33                           | 48                  |
| CPR-0630-143RX | 6                        | 30                     | 40                  | 40                           | 53                  | 43                           | 58                  |
| CPR-0718-143RX | 7                        | 18                     | 25                  | 23                           | 35                  | 26                           | 38                  |
| CPR-0725-143RX | 7                        | 25                     | 35                  | 30                           | 47                  | 36                           | 52                  |
| CPR-0730-143RX | 7                        | 30                     | 40                  | 40                           | 53                  | 44                           | 60                  |
| CPR-0820-143RX | 8                        | 20                     | 25                  | 25                           | 36                  | 27                           | 40                  |
| CPR-0825-143RX | 8                        | 25                     | 35                  | 30                           | 49                  | 38                           | 54                  |
| CPR-0830-143RX | 8                        | 30                     | 40                  | 40                           | 55                  | 45                           | 61                  |
| CPR-0840-143RX | 8                        | 40                     | 47                  | 50                           | 67                  | 60                           | 75                  |
| CPR-0920-143RX | 9                        | 20                     | 33                  | 30                           | 45                  | 33                           | 48                  |
| CPR-0930-143RX | 9                        | 30                     | 40                  | 40                           | 55                  | 45                           | 60                  |
| CPR-1020-143RX | 10                       | 20                     | 35                  | 30                           | 45                  | 35                           | 50                  |
| CPR-1030-143RX | 10                       | 30                     | 43                  | 40                           | 55                  | 45                           | 60                  |

### Features

Smallest stent pore size

5F Rx delivery system for all stent sizes

Self-tapering braided stent design

Stent can be recaptured after partial deployment

Closed cell stent with wall apposition and flexibility of an open cell design

### Benefits


Tacks down/contains plaque reducing embolic release

Very low profile

Optimize inventory with less model configurations

Improved accuracy of stent placement/landing zone

Stent conforms and adapts to all carotid anatomy



| <b>PHIL™</b><br>Non-adhesive Liquid Embolic System   |               |           |                     |
|--|---------------|-----------|---------------------|
| Product Code   | Concentration | Viscosity | Content Per Syringe |
| 1 pre-filled syringe of PHIL, 1 pre-filled syringe of DMSO and microcatheter hub adaptors* |               |           |                     |
| LEN10250   | 25%           | Low       | 1mL                 |
| LEN10300   | 30%           | Medium    | 1mL                 |
| LEN10350   | 35%           | High      | 1mL                 |

\*Available for Headway® DUO Microcatheter, Headway®17 Microcatheter, Scepter C® and Scepter XC® Balloon Occlusion Catheters and BALT® SONIC Microcatheter.

## Features

### 3 Concentrations Available

### Pre-loaded Syringes

### Iodine Bonded Radiopacifier

### High Volume of Embolic Precipitate

## Benefits

To address different clinical needs involving embolization of low flow/high flow lesions or a need for distal/proximal penetration

No preparation required, ready to use

Homogenous radiopacity  
No need to shake  
No added metal for minimal artifact

More embolic capacity with less DMSO agent

## Microcatheter Compatibility

PHIL™ must be used with DMSO compatible catheters such as Headway® DUO Microcatheter, Headway®17 Microcatheter or Scepter Occlusion Balloon.

| Product Name              | Product Codes                 | Description                       | Working Length (cm) | Dead Space (cc) | Dead Space with Adaptor (cc) | OD Prox./Dist. (French) | Tip Markers |
|---------------------------|-------------------------------|-----------------------------------|---------------------|-----------------|------------------------------|-------------------------|-------------|
| <b>Scepter C®</b>         | BC0410C<br>BC0415C<br>BC0420C | Compliant Occlusion Balloon       | 150                 | 0.44            | 0.23                         | 2.8/2.1                 | 3           |
| <b>Scepter XC®</b>        | BC0411XC                      | X-tra Compliant Occlusion Balloon | 150                 | 0.44            | 0.23                         | 2.8/2.1                 | 3           |
| <b>Headway® 17</b>        | MC172150S<br>MC172150SX       | Headway 17 Microcatheter          | 150                 | 0.44            | 0.26                         | 2.4/1.7                 | 2           |
| <b>Headway® Duo 156cm</b> | MC162156S                     | Headway Duo Microcatheter         | 156                 | 0.34            | 0.24                         | 2.1/1.6                 | 2           |
| <b>Headway® Duo 167cm</b> | MC162167S                     | Headway Duo Microcatheter         | 167                 | 0.35            | 0.25                         | 2.1/1.3                 | 1           |

# Ordering Information

To place an order, contact a MicroVention Customer Service Representative at:

## Ordering Information Customer Service

### MicroVention, Inc.

Attn: Customer Service  
1311 Valencia Avenue  
Tustin, CA 92780

### N. Am. Toll Free

(800) 990 8368

### Local

(714) 247 8000

### Customer Service Fax

(714) 247 8007

### Web

www.microvention.com

Please call customer service at 714 247 8000 to be connected with your local representative.

## MicroVention Locations Tustin Worldwide Headquarters and Europe

### MicroVention, Inc.

#### Worldwide Headquarters

1311 Valencia Avenue  
Tustin, CA 92780 USA  
PH +1.714.247.8000  
PH 1.800.990.8368

### MicroVention UK Limited

Suite 3, The Barracks Building  
10 Cliffords Fort, North Shields  
Tyne and Wear, NE30 1JE  
United Kingdom  
PH +44 (0) 191 258 6777  
F +44 (0) 191 258 5999

### MicroVention Europe, S.A.R.L.

30 bis, rue du Vieil Abreuvoir  
78100 Saint-Germain-en-Laye  
France  
PH +33 (1) 39 21 77 46  
F +33 (1) 39 21 16 01

### MicroVention Deutschland GmbH

Hildebrandtstr. 4 F  
D-40215 Düsseldorf  
Germany  
PH +49 211 210 798-0  
F +49 211 210 798-29

### Ordering Information:

Customer Service Hours are 7:00 am - 4:30 pm PST. Orders may be placed by phone, fax or mail. If mailing, please mark Attn: Customer Service. All prices and product specifications are subject to change at the discretion of MicroVention, Inc. MicroVention, Inc. warrants its products to be free from defect in materials and workmanship.

### Shipping Terms:

Orders will be shipped FOB shipping point using Federal Express Standard Overnight service, or comparable, unless otherwise requested. All freight charges are prepaid and added to invoice.

### Payment Terms:

Net 30 days.

### Returned Goods Policy:

Authorization is required for all returns.

1. Request Returned Goods Authorization (RGA) number from Customer Service.
2. Once the authorization number has been obtained, Customer Service will provide all the information necessary for a proper return.
3. Return product to:  
MicroVention, Inc.  
RGA # \_\_\_\_\_  
1311 Valencia Avenue  
Tustin, CA 92780

MicroVention, Inc. will accept the return of any standard catalog item under the following conditions:

- Merchandise processed or shipped in error by MicroVention, Inc. – *Full Credit*
- Product failure due to defect in materials or workmanship. – *Full Credit or Replacement at no charge.*